

100 bp DNA Ladder

100 to 1000 bp log scale

Ready-to-use DNA ladder

Cat.-No.	Volume	Amount
M-214	500 µl	100 lanes

For *in vitro* use only

Quality guaranteed for 12 months

Store at -20°C, short term storage (up to 4 weeks) at 4°C possible

100 bp DNA Ladder, ready-to-use

DNA fragment mix in 10 mM Tris-HCl (pH 7.5), 10 mM EDTA, 10% (w/v) glycerol, orange G and xylene cyanol FF

Description

100 bp DNA Ladder consists of 10 DNA fragments ranging from 100 bp to 1000 bp and is supplied in ready-to-use format containing blue tracking dyes. It allows sizing and concentration estimate of DNA fragments on agarose gels generated by PCR or restriction digest. The fragment mix yields a ladder of the following double-stranded DNA fragments: 100, 200, 300, 400, **500**, 600, 700, 800, 900, 1000 bp DNA log scale ladders are well suited for combination if small and large fragments should be analyzed on one gel.


Concentration

105 ng/µl

Recommended Load

5 µl DNA ladder per lane

100 bp DNA Ladder



2.4%
Agarose