

High Fidelity Hot Start Pol

Heat-activatable DNA polymerase for high specificity and accuracy

Thermus spezieis, recombinant, *E. coli*

	Cat. No.	Size	Conc.
	HFHS_10KU	10 kU	2.5 units/μl
	HFHS_100KU	100 kU	2.5 units/μl

For *in vitro* use only

Quality guaranteed for 12 months

Store at -20°C, avoid frequent thawing and freezing

Description

High Fidelity Hot Start Pol is based on a blend of Taq DNA polymerase and a proofreading enzyme specially designed for highly accurate and efficient amplification. The additional hot-start function provides improved specificity and sensitivity when amplifying low-copy-number targets in complex backgrounds or when prolonged room-temperature set up is required. The polymerase activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of nonspecifically annealed primers and primer-dimer formation at low temperatures during PCR setup. The enzyme shows excellent results with extremely long (up to 30 kb), GC-rich or other difficult templates. The enzyme blend includes a highly processive 5'→3' DNA polymerase and possesses a 5'→3' polymerization-dependent exonuclease replacement activity. Its inherent 3'→5' exonuclease proofreading activity results in a greatly increased fidelity of DNA synthesis compared to Taq polymerase.

Fidelity of the enzyme

High Fidelity Pol is characterized by a 4-fold higher fidelity compared to Taq polymerase.

$$ER_{\text{High Fidelity Pol}} = 3.4 \times 10^{-6}$$

The error rate (ER) of a PCR reaction is calculated using the equation $ER = MF / (bp \times d)$, where MF is the mutation frequency, bp is the number of base pairs of the fragment and d is the number of doublings ($2^d = \text{amount of product} / \text{amount of template}$).

Unit definition

One unit is defined as the amount of the enzyme required to catalyze the incorporation of 10 nmol of dNTP into an acid-insoluble form in 30 minutes at 74°C.

Recommended PCR assay

50 μl PCR assay		
5 μl	10x High fidelity buffer	green cap
200 μM	each dNTP	
0.2-0.5 μM	forward Primer	

High Fidelity Pol (red cap)

2.5 units/μl High Fidelity Hot Start Polymerase in storage buffer

10x High Fidelity Buffer (green cap)

High Fidelity Pol

Thermostable DNA polymerase for high accuracy

Thermus spezieis, recombinant, *E. coli*



0.2-0.5 µM	reverse Primer	
1-100 ng	Template DNA	
0.5 µl (1.25 units)	High Fidelity Pol	red cap
Fill up to 50 µl	PCR grade H ₂ O	

Please note that it is essential to add the polymerase as last component.

Recommended thermocycling conditions

Initial denaturation	95°C	2 min	1x
Denaturation	95°C	20 sec	20-30x
Annealing ¹⁾	50-68°C	30 sec	
Elongation ²⁾	72°C	1 min / kbp	
Final elongation	72°C	1 min / kbp	1x

- 1) The annealing temperature depends on the melting temperature of the primers used.
- 2) For amplification of fragments longer than 5 kb the elongation temperature should be set to 68°C.

For optimal specificity and amplification an individual optimization of the recommended parameters may be necessary for each new primer-template pair.