

High Fidelity Pol

Thermostable DNA polymerase for high accuracy

Thermus spezies, recombinant, *E. coli*

| | Cat. No. | Size | Conc. |
|--|----------|--------|--------------------|
| | HF_10KU | 10 kU | 2.5 units/ μ l |
| | HF_100KU | 100 kU | 2.5 units/ μ l |

For *in vitro* use only
Quality guaranteed for 12 months
Store at -20°C, avoid frequent thawing and freezing

Description

High Fidelity Pol is based on a blend of Taq DNA polymerase and a proofreading enzyme specially designed for highly accurate and efficient amplification. It shows excellent results with extremely long (up to 30 kb), GC-rich or other difficult templates. The enzyme blend includes a highly processive 5'→3' DNA polymerase and possesses a 5'→3' polymerization-dependent exonuclease replacement activity. Its inherent 3'→5' exonuclease proofreading activity results in a greatly increased fidelity of DNA synthesis compared to Taq polymerase. The enzyme is highly purified and free of bacterial DNA.

Fidelity of the enzyme

High Fidelity Pol is characterized by a 4-fold higher fidelity compared to Taq polymerase.

$$ER_{\text{High Fidelity Pol}} = 3.4 \times 10^{-6}$$

The error rate (ER) of a PCR reaction is calculated using the equation $ER = MF / (bp \times d)$, where MF is the mutation frequency, bp is the number of base pairs of the fragment and d is the number of doublings ($2^d = \text{amount of product} / \text{amount of template}$).

Unit definition

One unit is defined as the amount of the enzyme required to catalyze the incorporation of 10 nmol of dNTP into an acid-insoluble form in 30 minutes at 74°C.

Recommended PCR assay

| 50 μ l PCR assay | | |
|--------------------------|----------------------------|-----------|
| 5 μ l | 10x High fidelity buffer | green cap |
| 200 μ M | each dNTP | |
| 0.2-0.5 μ M | forward Primer | |
| 0.2-0.5 μ M | reverse Primer | |
| 1-100 ng | Template DNA | |
| 0.5 μ l (1.25 units) | High Fidelity Pol | red cap |
| Fill up to 50 μ l | PCR grade H ₂ O | |

Please note that it is essential to add the polymerase as last component.

High Fidelity Pol (red cap)

2.5 units/ μ l High Fidelity Polymerase in storage buffer

10x High Fidelity Buffer (green cap)

High Fidelity Pol

Thermostable DNA polymerase for high accuracy

Thermus spezies, recombinant, *E. coli*



Recommended thermocycling conditions

| | | | |
|----------------------------|---------|-------------|--------|
| Initial denaturation | 95°C | 2 min | 1x |
| Denaturation | 95°C | 20 sec | 20-30x |
| Annealing ¹⁾ | 50-68°C | 30 sec | |
| Elongation ^{2,3)} | 72°C | 1 min / kbp | |
| Final elongation | 72°C | 1 min / kbp | 1x |

- 1) The annealing temperature depends on the melting temperature of the primers used.
- 2) For amplification of fragments longer than 5 kb the elongation temperature should be set to 68°C.
- 3) The elongation time depends on the length of the fragments to be amplified. A time of 1 min/kbp is recommended.

For optimal specificity and amplification an individual optimization of the recommended parameters may be necessary for each new primer-template pair.